




## VOC METHOD UPDATE SKC APPENDICES TO EPA METHOD TO-17

Year 2015  
www.skcinc.com

*This method update has been written by SKC as a guideline for users. The sampling apparatus specified in this SKC update reflects new technology that may not have been available at the time of the original publication. This method update by SKC has not been officially endorsed or approved by U.S. EPA.*

### DETERMINATION OF VOLATILE ORGANIC COMPOUNDS (VOCs) IN AMBIENT OR INDOOR AIR

- **SKC APPENDIX A**  
**Active Sampling Using SKC Solid Adsorbent Tubes**



- **SKC APPENDIX B**  
**Passive Sampling Using SKC ULTRA Series Diffusive Samplers**



## Scope

1.1 This document is written by SKC as appendices to U.S. EPA Compendium Method TO-17, *Determination of Volatile Organic Compounds in Ambient Air Using Active Sampling onto Sorbent Tubes* (published in 1997). **The samplers described, however, are suitable for sampling VOCs in ambient or indoor air.** As such, this document also provides an update of EPA Compendium Method IP-1B, *Determination of Volatile Organic Compounds in Indoor Air Using Solid Adsorbent Tubes* (published in 1989).

1.2 Method TO-17 specifies the use of a solid adsorbent tube for sampling followed by thermal desorption (TD) and gas chromatography (GC) analysis. **SKC Appendix A**, presented in this document, provides a listing of solid sorbent tubes available from SKC for the compounds referenced in the original TO-17 method.

1.3 **SKC Appendix B** describes a passive sampling method for VOCs using diffusive samplers manufactured by SKC. With these samplers, VOCs diffuse onto the sorbent at a fixed, defined rate determined by the specific chemical and the geometry of the samplers. The sorbent is transferred from the sampler into a tube for thermal desorption and gas chromatographic analysis.

1.4 Both the active and passive sampling methods are capable of measuring target compounds in the low parts-per-billion (ppb) and parts-per-trillion (ppt) concentration range.

## 2. Significance

2.1 Accurate risk assessment of human and ecological exposure to toxic VOCs is an important goal of the U.S. EPA. Of the 189 hazardous air pollutants (HAPs) listed in Title III of the Clean Air Act Amendments (CAAA) of 1990, 97 are VOCs. Similarly, in the indoor environment, VOCs are high priority pollutants. Construction of more tightly sealed buildings and the use of synthetic building materials and furnishings can significantly increase the indoor levels of many VOCs. In addition, VOCs can migrate from buried waste and/or contaminated groundwater through subsurface solid into the air space of overlying buildings (vapor intrusion). These vapors may then accumulate, thereby reducing the indoor air quality and posing a health risk to building occupants. VOC vapor intrusion is a hot-button issue in the field of environmental science and regulatory policy. Therefore, measuring VOC levels is an essential part of evaluating potential health threats and identifying appropriate controls.

2.2 Historically, stainless steel canisters have been used to measure VOC exposures in ambient air. Canisters, however, do not have the size and portability for micro-environmental and human subject studies. In addition, the cost and logistics of canister sampling can be quite burdensome. [\(8\)](#)

**2.3** Solid sorbent samplers offer users a simpler option for measuring VOCs. EPA Compendium Method TO-17 describes several advantages of solid sorbent samplers:

- The small size and light weight of the sorbent sampler and attendant equipment
- The availability of a large selection of sorbents to match the target compounds
- The commercial availability of thermal desorption systems to release the sample from the sorbent and into the analytical system
- The placement of the sorbent packing as the first element in the sampling train, thus reducing the possibility of contamination or losses from upstream elements
- The large amount of literature on the use of sorbent sampling and thermal desorption for monitoring workplace air, particularly from the Health and Safety Executive (HSE) in the United Kingdom (U.K.).

**2.4** The procedures described in these appendices provide the user with a choice of methodologies for the sampling and analysis of VOCs in ambient or indoor air using solid sorbent technologies: active sampling with a pump/sorbent tube or passive sampling with a diffusive sorbent sampler.

# DETERMINATION OF VOLATILE ORGANIC COMPOUNDS IN AMBIENT AIR USING ACTIVE SAMPLING ONTO SORBENT TUBES

## EPA METHOD TO-17 HIGHLIGHTS

**Sampler** - Glass or stainless steel tubes of various lengths and outer diameters (OD) with the central portion packed with  $\geq 200$  mg of solid adsorbent material(s). Tubes should be individually numbered and show the direction of flow.

**Flow rates** - Two tubes<sup>§</sup> with independent flow control should be run in parallel to achieve sample volumes of 1 and 4 liters<sup>‡</sup>. Typical flow rates include:

- 16.7 ml/min to collect 1 liter of air in 1 hour
- 66.7 ml/min to collect 4 liters of air in 1 hour

§ *EPA requires the use of distributed volume pairs for monitoring to ensure high quality data. However, in situations where acceptable data have been routinely obtained through use of distributed volume pairs and the ambient air is considered well characterized, cost considerations may warrant single-tube sampling. Any attendant risk to data quality objectives is the responsibility of the individual.*

‡ *Appropriate proportional scaling of these volumes to fit the target list and monitoring objectives is acceptable.*

**Sample time** - 1 to 3 hours

**Sample temperature** - The temperature of the sorbent tube must be the same (and not lower than) ambient temperature when sampling or moisture will be retained via condensation.

**Sample storage** - Use sample tubes within 30 days after conditioning with heat and an inert gas. After sampling, place the tubes in a clean, cool ( $< 4$  C), and organic solvent-free environment. Generally, analysis should occur within 30 days of sample collection as artifacts may be generated during sample storage. Some compounds, such as limonene and labile sulfur or nitrogen-containing volatiles, will require analysis within one week. Tubes with multiple sorbent beds should also be analyzed as soon as possible to avoid sample recovery errors. Blank tubes should be stored in a manner similar to the sample tubes to check background levels of the sorbent.

**Analysis** - Thermal desorption followed by gas chromatography

**Background levels of samplers** - Tubes packed with sorbent should be conditioned using elevated temperatures and a flow of ultra-pure inert gas. Table 2 in [EPA Method TO-17](#) provides guidelines on the recommended temperature and gas flow for conditioning. Analyze a conditioned tube to obtain a blank profile for target compounds.

**Capacity of samplers** - Method TO-17, Appendix 1, provides a list of safe sample volumes (SSV) for target compounds.

**Typical limits of detection** - A method detection limit of < 0.5 ppb is a required performance criterion for Method TO-17. Detection limits will vary depending on several factors.

**Precision** - Duplicate (analytical) precision should be within 20%.

**Audit accuracy** - The degree of agreement with a known audit standard that should be within 30%

**Distributed volume pairs** - The absolute value of the relative difference between the two tube samples should be within 25%.

VOC METHOD UPDATE  
**DETERMINATION OF VOLATILE ORGANIC COMPOUNDS IN  
 AMBIENT OR INDOOR AIR**

**EPA METHOD TO-17 – SKC APPENDIX A  
 (Active Solid Sorbent Tube)**

EPA Compendium Method TO-17 can be found in its entirety online at [www.epa.gov/ttnamti1/files/ambient/airtox/to-17r.pdf](http://www.epa.gov/ttnamti1/files/ambient/airtox/to-17r.pdf). SKC Appendix A lists SKC sorbent tubes that are suitable for target compounds, type of sorbent contained in each tube, and the tube material. For further details on SKC thermal desorption tubes, visit [www.skcinc.com/prod/226-356.asp](http://www.skcinc.com/prod/226-356.asp) or see the SKC Catalog and Sampling Guide.

**HYDROCARBONS**

Aliphatic, aromatic, and cyclic hydrocarbons less volatile than ethane and more volatile than nC<sub>20</sub>

<b>Target Compound</b>	<b>Sorbent Type*</b>			<b>Tube</b> (G=Glass) (SS=Stainless Steel)	<b>SKC Cat. No.</b>
	<b>Bed 1</b>	<b>Bed 2</b>	<b>Bed 3</b>		
n-Butane n-Pentane	Anasorb® GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
n-Pentane	Anasorb GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
	Chromosorb® 106	-----	-----	SS	226-358
n-Hexane	Tenax® GR	Anasorb GCB1	-----	G	226-345
	Anasorb GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Tenax GR	Anasorb GCB1	-----	SS	226-348
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
	Anasorb GCB1	-----	-----	SS	226-356
Chromosorb 106	-----	-----	SS	226-358	

\* Anasorb GCB1 is equivalent to Carbopack® B; Anasorb GCB2 is equivalent to Carbopack C.

## HYDROCARBONS (continued)

<b>Target Compound</b>	<b>Sorbent Type*</b>			<b>Tube</b> (G=Glass) (SS=Stainless Steel)	<b>SKC Cat. No.</b>
	<b>Bed 1</b>	<b>Bed 2</b>	<b>Bed 3</b>		
Benzene n-Heptane Toluene n-Octane Ethylbenzene Xylenes n-Nonane 1,3,5-Trimethyl benzene n-Decane	Tenax GR	Anasorb GCB1	-----	G	226-345
	Anasorb GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Tenax GR	Anasorb GCB1	-----	SS	226-348
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
	Anasorb GCB1	-----	-----	SS	226-356
	Tenax TA	-----	-----	SS	226-357
	Chromosorb 106	-----	-----	SS	226-358
	Tenax TA	-----	-----	G	226-360
Styrene Isopropylbenzene n-Propylbenzene 1-Methyl-3- ethylbenzene 1-Methyl-4- ethylbenzene Methyl styrene Methyl-2-ethyl benzene 1,2,4-Trimethyl benzene 1,2,3-Trimethyl benzene n-Undecane	Tenax GR	Anasorb GCB1	-----	G	226-345
	Anasorb GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Tenax GR	Anasorb GCB1	-----	SS	226-348
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
	Anasorb GCB1	-----	-----	SS	226-356
	Tenax TA	-----	-----	SS	226-357
	Tenax TA	-----	-----	G	226-360
	n-Dodecane	Tenax GR	Anasorb GCB1	-----	G
Anasorb GCB2		Anasorb GCB1	Carbosieve S-III	G	226-347
Tenax GR		Anasorb GCB1	-----	SS	226-348
Anasorb GCB2		Anasorb GCB1	Carbosieve S-III	SS	226-350
Anasorb GCB1		-----	-----	SS	226-356
Tenax TA		-----	-----	SS	226-357
Tenax TA		-----	-----	G	226-360

\* Anasorb GCB1 is equivalent to Carbo-pack B; Anasorb GCB2 is equivalent to Carbo-pack C.

## HALOGENATED HYDROCARBONS INCLUDING PCBs

Aliphatic, aromatic, and cyclic halogenated hydrocarbons more volatile than nC<sub>20</sub>

Target Compound	Sorbent Type*			Tube (G=Glass) (SS=Stainless Steel)	SKC Cat. No.
	Bed 1	Bed 2	Bed 3		
Dichloromethane	Anasorb GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
1,2-Dichloroethane	Tenax GR	Anasorb GCB1	-----	G	226-345
	Anasorb GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Tenax GR	Anasorb GCB1	-----	SS	226-348
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
	Anasorb GCB1	-----	-----	SS	226-356
	Tenax TA	-----	-----	SS	226-357
	Chromosorb 106	-----	-----	SS	226-358
Tenax TA	-----	-----	G	226-360	
1,1,1-Trichloroethane	Tenax GR	Anasorb GCB1	-----	G	226-345
	Anasorb GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Tenax GR	Anasorb GCB1	-----	SS	226-348
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
	Anasorb GCB1	-----	-----	SS	226-356
	Chromosorb 106	-----	-----	SS	226-358
Carbon Tetrachloride Trichloroethylene 1,1,2-Trichloroethane Tetrachloroethylene Chlorobenzene 1,1,1,2-Tetrachloroethane 1,1,2,2-Tetrachloroethane	Tenax GR	Anasorb GCB1	-----	G	226-345
	Anasorb GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Tenax GR	Anasorb GCB1	-----	SS	226-348
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
	Anasorb GCB1	-----	-----	SS	226-356
	Tenax TA	-----	-----	SS	226-357
	Chromosorb 106	-----	-----	SS	226-358
	Tenax TA	-----	-----	G	226-360

\* Anasorb GCB1 is equivalent to Carbo-pack B; Anasorb GCB2 is equivalent to Carbo-pack C.

## ALCOHOLS

Alcohols more volatile than nC<sub>20</sub> and sufficiently stable to be analyzed by conventional GC techniques

Target Compound	Sorbent Type*			Tube (G=Glass) (SS=Stainless Steel)	SKC Cat. No.
	Bed 1	Bed 2	Bed 3		
Methanol Ethanol	Anasorb GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
n-Propanol Isopropanol	Tenax GR	Anasorb GCB1	-----	G	226-345
	Anasorb GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Tenax GR	Anasorb GCB1	-----	SS	226-348
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
	Anasorb GCB1	-----	-----	SS	226-356
Chromosorb 106	-----	-----	SS	226-358	
n-Butanol iso-Butanol	Tenax GR	Anasorb GCB1	-----	G	226-345
	Anasorb GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Tenax GR	Anasorb GCB1	-----	SS	226-348
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
	Anasorb GCB1	-----	-----	SS	226-356
	Tenax TA	-----	-----	SS	226-357
	Chromosorb 106	-----	-----	SS	226-358
Tenax TA	-----	-----	G	226-360	
Octanol	Tenax GR	Anasorb GCB1	-----	G	226-345
	Anasorb GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Tenax GR	Anasorb GCB1	-----	SS	226-348
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
	Tenax TA	-----	-----	SS	226-357
	Tenax TA	-----	-----	G	226-360

\* Anasorb GCB1 is equivalent to Carboxpack B; Anasorb GCB2 is equivalent to Carboxpack C.

## ESTERS AND GLYCOL ETHERS

Esters and glycol ethers more volatile than nC<sub>20</sub> and sufficiently stable to be analyzed by conventional GC techniques

<b>Target Compound</b>	<b>Sorbent Type*</b>			<b>Tube</b> (G=Glass) (SS=Stainless Steel)	<b>SKC Cat. No.</b>
	<b>Bed 1</b>	<b>Bed 2</b>	<b>Bed 3</b>		
Methyl acetate Methyl-t-butyl ether t-Butyl acetate	Tenax GR	Anasorb GCB1	-----	G	226-345
	Anasorb GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Tenax GR	Anasorb GCB1	-----	SS	226-348
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
	Anasorb GCB1	-----	-----	SS	226-356
	Chromosorb 106	-----	-----	SS	226-358
Ethyl acetate Propyl acetate Isopropyl acetate Butyl acetate Isobutyl acetate Methyl acrylate Ethyl acrylate Methyl metha- crylate Methoxyethanol Ethoxyethanol Butoxyethanol Methoxypropanol Methoxy- ethylacetate Ethoxyethylacetate Butoxyethylacetate	Tenax GR	Anasorb GCB1	-----	G	226-345
	Anasorb GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Tenax GR	Anasorb GCB1	-----	SS	226-348
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
	Anasorb GCB1	-----	-----	SS	226-356
	Tenax TA	-----	-----	SS	226-357
	Chromosorb 106	-----	-----	SS	226-358
	Tenax TA	-----	-----	G	226-360

\* Anasorb GCB1 is equivalent to Carbopack B; Anasorb GCB2 is equivalent to Carbopack C.

## ALDEHYDES AND KETONES

Aldehydes and ketones more volatile than nC<sub>20</sub> and sufficiently stable to be analyzed by conventional GC techniques

<b>Target Compound</b>	<b>Sorbent Type*</b>			<b>Tube</b> (G=Glass) (SS=Stainless Steel)	<b>SKC Cat. No.</b>
	<b>Bed 1</b>	<b>Bed 2</b>	<b>Bed 3</b>		
Acetone	Anasorb GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
	Chromosorb 106	-----	-----	SS	226-358
2-Butanone Methylisobutyl- ketone Cyclohexanone	Tenax GR	Anasorb GCB1	-----	G	226-345
	Anasorb GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Tenax GR	Anasorb GCB1	-----	SS	226-348
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
	Anasorb GCB1	-----	-----	SS	226-356
	Tenax TA	-----	-----	SS	226-357
	Chromosorb 106	-----	-----	SS	226-358
Tenax TA	-----	-----	G	226-360	
n-Butanal	Tenax GR	Anasorb GCB1	-----	G	226-345
	Anasorb GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Tenax GR	Anasorb GCB1	-----	SS	226-348
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
	Anasorb GCB1	-----	-----	SS	226-356
	Chromosorb 106	-----	-----	SS	226-358
3,5,5-Trimethyl- cyclohex-2-enone	Tenax GR	Anasorb GCB1	-----	G	226-345
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Tenax GR	Anasorb GCB1	-----	SS	226-348
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
	Anasorb GCB1	-----	-----	SS	226-356
	Tenax TA	-----	-----	SS	226-357

\* Anasorb GCB1 is equivalent to Carbowack B; Anasorb GCB2 is equivalent to Carbowack C.

## ALDEHYDES AND KETONES (continued)

Target Compound	Sorbent Type*			Tube (G=Glass) (SS=Stainless Steel)	SKC Cat. No.
	Bed 1	Bed 2	Bed 3		
Furfural	Tenax GR	Anasorb GCB1	-----	G	226-345
	Anasorb GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Tenax GR	Anasorb GCB1	-----	SS	226-348
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
	Anasorb GCB1	-----	-----	SS	226-356
	Tenax TA	-----	-----	SS	226-357
Tenax TA	-----	-----	G	226-360	

\* Anasorb GCB1 is equivalent to Carbopack B; Anasorb GCB2 is equivalent to Carbopack C.

## MISCELLANEOUS VOCs

These are most VOCs in air. Generally compatible with all organics less volatile than ethane, more volatile than nC<sub>20</sub>, and sufficiently stable to be analyzed by conventional GC techniques.

Target Compound	Sorbent Type*			Tube (G=Glass) (SS=Stainless Steel)	SKC Cat. No.
	Bed 1	Bed 2	Bed 3		
Acetonitrile	Anasorb GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
Acrylonitrile Propionitrile Acetic acid	Tenax GR	Anasorb GCB1	-----	G	226-345
	Anasorb GCB1	Carbosieve S-III		G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Tenax GR	Anasorb GCB1	-----	SS	226-348
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
Anasorb GCB1	-----	-----	SS	226-356	

\* Anasorb GCB1 is equivalent to Carbopack B; Anasorb GCB2 is equivalent to Carbopack C.

### MISCELLANEOUS VOCs (continued)

<b>Target Compound</b>	<b>Sorbent Type*</b>			<b>Tube</b> (G=Glass) (SS=Stainless Steel)	<b>SKC Cat. No.</b>
	<b>Bed 1</b>	<b>Bed 2</b>	<b>Bed 3</b>		
Maleic anhydride Pyridine Aniline	Tenax GR	Anasorb GCB1	-----	G	226-345
	Anasorb GCB1	Carbosieve S-III	-----	G	226-346
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Tenax GR	Anasorb GCB1	-----	SS	226-348
	Anasorb GCB1	Carbosieve S-III	-----	SS	226-349
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	SS	226-350
	Anasorb GCB1	-----	-----	SS	226-356
	Tenax TA	-----	-----	SS	226-357
	Chromosorb 106	-----	-----	SS	226-358
	Tenax TA	-----	-----	G	226-360
Nitrobenzene	Tenax GR	Anasorb GCB1	-----	G	226-345
	Anasorb GCB2	Anasorb GCB1	Carbosieve S-III	G	226-347
	Tenax GR	Anasorb GCB1	-----	SS	226-348
	Anasorb GCB1	Anasorb GCB1	Anasorb GCB1	SS	226-350
	Tenax TA	-----	-----	SS	226-357
	Tenax TA	-----	-----	G	226-360
Phenol	Tenax GR	Anasorb GCB1	-----	G	226-345
	Tenax GR	Anasorb GCB1	-----	SS	226-348
	Tenax TA	-----	-----	SS	226-357
	Tenax TA	-----	-----	G	226-360

\* Anasorb GCB1 is equivalent to Carbopack B; Anasorb GCB2 is equivalent to Carbopack C.

**VOC METHOD UPDATE  
DETERMINATION OF VOLATILE ORGANIC COMPOUNDS IN  
AMBIENT AND INDOOR AIR**

**EPA METHOD TO-17 – SKC APPENDIX B  
(Passive [Diffusive] Sampler)**

## *SKC APPENDIX B HIGHLIGHTS*

**Sampler** - A nylon sampler housing with a number of inlet holes is loaded with one of several sorbents. Sampling is started by removing or sliding the cover to expose the holes to the air and is stopped by replacing the cover (SKC Inc, Eighty Four, PA; ULTRA Series Passive Samplers, Cat. No. 690 Series). With these samplers, VOCs diffuse onto the sorbent at a fixed, defined rate determined by the specific chemical and the geometry of the samplers.

**Sampling/Uptake rate** - Varies by compound. SKC and the U.S. OSHA Testing Lab have experimentally verified the sampling or uptake rates of many chemicals using SKC passive samplers. See the online [Passive Sampling Guide](#) for compound-specific sampling/uptake rates. Despite differences in the outward appearance of samplers, the sampler housings maintain the same geometry; therefore, the sampling rates are equivalent for all models provided similar sorbents are used. Additional information on sampling rates using Chromosorb 106 is available. [\(5\)](#)

**Face velocity and orientation** - The SKC ULTRA Series passive samplers will work in face or wind velocities of 5 to 200 cm/sec. Wind velocities < 5 cm/sec may result in reduced sampling rates. Operation is proficient in a parallel or perpendicular orientation to wind. For more information, go to [www.skcinc.com](http://www.skcinc.com) and enter "ULTRA" in the Google search box.

**Sample time** – 8 hours to 30 days

**Sampling conditions** - Temperature range: 10 to 40 C and relative humidity range: 20 to 80% RH (25 C)

**Sample storage - Before use:** Store at ≤ 39.2 F (4 C) for up to 30 days to reduce the formation of artifacts and spurious peaks in the sorbent. **After use:** Store at ≤ 39.2 F (4 C) and analyze within 21 days. Store in a clean and organic solvent-free environment.

**Analysis** - Sorbent is transferred directly from the sampler housing to a thermal desorption tube (0.25-inch OD and minimum 3.5-inch length). The tube is sealed until analysis. The VOCs thermally desorbed from the sorbent are analyzed by gas chromatography, typically with mass spectrometry (MS).

**Background levels of samplers** - Typically less than 25 ng per passive sampler

# DETERMINATION OF VOLATILE ORGANIC COMPOUNDS IN AMBIENT AND INDOOR AIR

## EPA METHOD TO-17 – SKC APPENDIX B (Passive [Diffusive] Sampler)

### TABLE OF CONTENTS

	<u>Page</u>
<a href="#">1.</a> Scope .....	16
<a href="#">2.</a> Applicable Documents .....	19
<a href="#">3.</a> Summary of Method .....	19
<a href="#">4.</a> Significance .....	20
<a href="#">5.</a> Definitions.....	20
<a href="#">6.</a> Interferences .....	20
<a href="#">7.</a> Apparatus .....	20
<a href="#">8.</a> Reagents and Materials.....	21
<a href="#">9.</a> Preparation of Reagents.....	21
<a href="#">10.</a> Sample Collection.....	21
<a href="#">11.</a> Analytical Procedure.....	21
<a href="#">12.</a> Method Detection Limits .....	22
<a href="#">13.</a> References .....	22
<a href="#">Figures</a> .....	24

## 1. Scope

**1.1** This method describes a procedure for the determination of VOCs in ambient or indoor air in the low-ppb and ppt range. The procedure has been found to be suitable for 8-hour to 30-day indoor/ambient air studies targeting medium-to-high boiling compounds.

**1.2** This method uses the same thermal desorption techniques and gas chromatographic analysis specified in active sampling Method TO-17. In this case, however, VOCs are collected through a passive process of controlled diffusion through inlet holes of a sampler loaded with solid sorbent (SKC Inc., Eighty Four, PA; ULTRA Series Passive Samplers, Cat. No. 690 Series). Samplers are available with five different sorbents. At the laboratory, sorbent is transferred to a thermal desorption tube with an outer diameter of 0.25 inch and a minimum length of 3.5 inches (see *Operating Instructions at [www.skcinc.com/instructions/40111.pdf](http://www.skcinc.com/instructions/40111.pdf)*). The tube is sealed until analysis. VOCs are thermally desorbed from the sorbent and analyzed by gas chromatography, typically with mass spectrometry.

**1.3** The uptake rate varies by compound. Uptake rates and other critical sampling parameters have been experimentally verified for many VOCs by SKC and other agencies. **(1-2)** See the [online Passive Sampling Guide](#) for compound-specific sampling/uptake rates. Despite differences in the outward appearance of samplers, the sampler housings maintain the same geometry; therefore, the sampling rates are equivalent for all models provided similar sorbents are used. Only reverse diffusion may differ based on the sorbent and the chemical being sampled.

**1.3.1** SKC ULTRA Series passive samplers will work in face or wind velocities of 5 to 200 cm/sec. Operation is proficient in a parallel or perpendicular orientation to wind. When face velocities are < 5 cm/sec, the sampling rate may be reduced (commonly called “zero velocity sampling rates”). These lower velocity rates are typically found when sampling inside homes or buildings. For more information, go to [www.skcinc.com](http://www.skcinc.com) and enter “ULTRA” in the Google search box.

**1.4** Researchers at the Sahlgrenska Academy at Goteborg University and the National Institute of Working Life in Sweden have also validated the use of SKC passive sampler for 1,3-butadiene and benzene in ambient air for eight hours up to seven days. (3)

**1.5** The passive samplers are marketed as the ULTRA Series Passive Samplers for thermal desorption and are manufactured by SKC Inc. (Eighty Four, PA; ULTRA Series Passive Samplers, Cat. No. 690 Series). Up to five sorbents are available.

**1.6** This method may involve hazardous materials, operations, and equipment. This method does not purport to address all the safety problems associated with its use. It is the user's responsibility to develop and implement appropriate safety and health practices and to determine the applicability of regulatory limitations prior to use.

## **2. Applicable Documents**

### **2.1 ASTM International Standards**

*D1356 Definition of Terms Relating to Atmospheric Sampling and Analysis*

*E682 Practice for Liquid Chromatography Terms and Relationships*

### **2.2 Other Documents**

- *Compendium of Methods for the Determination of Air Pollutants in Indoor Air*
- *Compendium of Methods for the Determination of Organic Compounds in Ambient Air*
- *Methods for the Determination of Hazardous Substances by the Health and Safety Laboratory - MDHS 88 Volatile Organic Compounds in Air*
- *Performance of SKC ULTRA Passive Samplers by the U.S. Department of Labor*

## **3. Summary of Method**

**3.1** Sorbent is loaded in a sampler housing that contains multiple inlet holes. VOC vapor diffuses into the sampler at an experimentally determined rate and is collected onto the solid sorbent material. The samplers have been found to be suitable for 8 hours up to 30 days of indoor/ambient air sampling targeting medium-to-high boiling compounds.

**3.2** After sampling, samplers are closed with the sampler covers (see *Operating Instructions* [at www.skcinc.com/instructions/40111.pdf](http://www.skcinc.com/instructions/40111.pdf)). Samplers are sealed in the provided pouches and packed for shipment to a laboratory. Store samples at  $\leq 39.2$  F (4 C) in a clean and organic solvent-free environment followed by analysis within 21 days.

**3.3** In the analytical laboratory, the sorbent is transferred from the sampler housing directly into a thermal desorption tube. The tube is then sealed with PTFE caps or Swagelok® fittings. Thermal desorption tubes must be an outer diameter of 0.25 inch and a minimum length of 3.5 inches, such as those manufactured by Perkin-Elmer.

**3.4** The VOCs adsorbed onto the sorbent material are thermally desorbed and analyzed by GC or GC/MS. The blank samples are likewise desorbed and analyzed.

**3.5** VOCs in the sample are identified and quantified by comparison with calibration curves for the analytes of interest.

## **4. Significance**

**4.1** This method uses a passive sampling system with an uptake rate that depends on the diffusion coefficient of the VOC of interest, and that is proportional to the cross-sectional area of the opening of the sampler and inversely proportional to the length of the diffusion zone of the sampler. SKC and other agencies have experimentally verified sampling/uptake rates for various compounds. The samplers are available commercially from SKC Inc. as the ULTRA Series Passive Samplers (see [www.skcinc.com/prod/690-101.asp](http://www.skcinc.com/prod/690-101.asp) for details). The sorbent material must be used within 30 days or repurged by the user to ensure acceptable background levels. Typical background levels of VOCs **are less than 25 ng** per sampler. The portable samplers allow for employment of this method in general areas with some minimal air movement or in micro-environments.

**4.2** Thermal desorption and subsequent GC or GC/MS analysis provide a very accurate measure of low-level VOC concentrations.

## 5. Definitions

**Note:** Definitions used in this document and any user-prepared SOPs should be consistent with ASTM International Methods D1356 and E682. All pertinent abbreviations and symbols are defined within this document at point of use.

## 6. Interferences

**6.1** The ULTRA Series Passive Samplers allow the user to select the best sorbent for the target compounds and transfer sorbent easily for analysis. The Marine Corps selected Tenax TA for the measurement of chemical warfare agents. (4) U.S. OSHA reported on the use of Carboxen 1016, Carbotrap Z<sup>®</sup>, and Chromosorb 106 for a mixture containing 20 solvent analytes. (5) SKC has conducted studies on Tenax TA and Anasorb GCB1 sorbents (6) and Carbograph 5 TD sorbent. (7)

**6.2** Interferences are not a problem if the sorbent selection and GC parameters are properly managed.

## 7. Apparatus

### 7.1 Sampling

#### 7.1.1 Passive Diffusive Samplers

**ULTRA** (see *Figure 1*) – A housing loaded with sorbent with inlet holes (SKC Inc., Eighty Four, PA; Cat. Nos. 690-101, -103, -104, -105, and -106) that allows for controlled diffusion and collection of VOCs. The sampler contains discrete compartments for sample sorbent and blank/correction sorbent. A cover slides for opening and closing the sampler. Sample sorbent is transferred to a thermal desorption tube.

Uptake rate is dependent on the sorbent, the target compound, and sampling conditions. The housing has a clip for positioning during sampling. SKC manufactures ULTRA Series samplers to meet the following criteria:

- Total VOC concentration: < 0.025 µg (25 ng) per sampler or vial

**7.1.2** The samplers are shipped from the manufacturer in a pouch, which can be used after sampling for shipping the samplers or sorbent vials in an appropriate mailer to a laboratory.

**7.1.3** Sorbent is typically conditioned by packing it inside a tube using elevated temperatures and a flow of ultra-pure inert gas. EPA Method TO-17, Table 2 ([www.epa.gov/ttnamti1/files/ambient/airtox/to-17r.pdf](http://www.epa.gov/ttnamti1/files/ambient/airtox/to-17r.pdf)) provides guidelines on the recommended temperature and gas flow for conditioning.

**7.1.4** Use sorbent within 30 days after conditioning with heat and an inert gas.

## **7.2 Analysis**

**7.2.1** After sampling, place sampler with sorbent or sorbent vial in a clean, cool ( $\leq 4$  C), and organic solvent-free environment. Generally, analysis should occur within 21 days of sample collection as artifacts may be generated during sample storage. Some compounds, such as limonene and labile sulfur or nitrogen-containing volatiles, will require analysis within one week. Blank samplers or sorbent vials (from the same lot as sample sorbent) should be stored in a manner similar to the samples to check background levels of the sorbent.

**7.2.2** Analytes should be desorbed from the sorbent by back flushing the tubes with heat and an inert gas (typically 150 to 400 C for 5 to 15 minutes with a carrier gas flow of 30 to 100 ml/min). The VOCs thermally desorbed from the sorbent are analyzed by gas chromatography, typically with mass spectrometry. *See sample chromatogram in Figure 2.*

## **8. Reagents and Materials**

### **9. Preparation of Reagents**

Complete information for Sections 8 and 9 can be obtained in EPA Method TO-17 at [www.epa.gov/ttn/amtic/airtox.html](http://www.epa.gov/ttn/amtic/airtox.html).

## **10. Sample Collection**

**10.1** See *Operating Instructions* at [www.skinc.com/instructions/40111.pdf](http://www.skinc.com/instructions/40111.pdf).

## **11. Analytical Procedure**

### **11.1 Sample Preparation and Desorption**

**11.1.1** See Analysis Instructions at [www.skinc.com/](http://www.skinc.com/); search on 40125 or 40127.

**11.1.2** Analytes should be desorbed from the sorbent by back flushing the tubes with heat and an inert gas (typically 150 to 400 C for 5 to 15 minutes with a carrier gas flow of 30 to 100 ml/min). The VOCs thermally desorbed from the sorbent are analyzed by gas chromatography with the detector recommended for the compound(s) of interest, typically mass spectrometry.

See EPA Method TO-17 at [www.epa.gov/ttn/amtic/airtox.html](http://www.epa.gov/ttn/amtic/airtox.html) for complete information on analysis, calibration, and quality assurance procedures.


## 12. Method Detection Limits

12.1 SKC researchers used the following procedure to determine detection limits.


12.1.1 Reporting limits were determined by running seven replicates of a standard at a level near the detection limit of the analyte. The standard deviation of the mean of these values was calculated and if the relative standard deviation was less than 25%, this was stated as the method reporting limit.

## 13. References

1. **SKC Inc.:** *Passive (Diffusive) Sampler Research Reports*. Eighty Four, PA. See online at [www.skcinc.com/reports.asp](http://www.skcinc.com/reports.asp).
2. **U.S. Department of Labor:** *OSHA Analytical Methods Manual*. Salt Lake City, Utah: U.S. Department of Labor, Occupational Safety and Health Division Directorate for Technical Support, 1990. See online at [www.osha.gov/](http://www.osha.gov/)
3. **Strandberg, B., Olsson, K., Levin, J., et. al.:** "[Evaluation of Two Types of Diffusive Samplers and Adsorbents for Measuring 1,3-Butadiene and Benzene in Air](#)," *Atmospheric Environment*, Vol. 39, Issue 22, July 2005, pp. 4101-4110
4. **Hendricks, W.D., Roberts, J.L., Schultz, G.R.:** "Feasibility of Diffusive Sampling to Monitor U.S. Military Personnel for Exposure to Toxic Chemical Substances," Salt Lake City, Utah: U.S. Department of Labor, Occupational Safety and Health Division Directorate for Technical Support, 2001
5. **U.S. Department of Labor:** "Performance of SKC ULTRA Passive Samplers Containing Carboxen 1016, Carbotrap Z, or Chromosorb 106 When Challenged with a Mixture Containing Twenty of OSHA SLTC's Top Solvent Analytes," Salt Lake City, Utah: U.S. Department of Labor, Occupational Safety and Health Division Directorate for Technical Support, 2003. See online at [www.osha.gov/dts/sltc/methods/studies/chrom\\_106\\_ultra/chrom\\_106\\_ultra.html](http://www.osha.gov/dts/sltc/methods/studies/chrom_106_ultra/chrom_106_ultra.html).
6. **Kuhlman, C., Coyne, L.:** "Feasibility of a New Diffusive Sampler to Measure Sub-ppb Levels of VOCs," SKC Inc., presented at the AWMA Symposium on Air Quality Measurement Methods and Technology, April 2005
7. **Coyne, L., Kuhlman, C., Chada, J.:** "Validation of the 575-001 and the ULTRA Series Diffusive Samplers: Long-term Sampling in Indoor and Ambient Air Environments," [www.skcinc.com/pdf/1812.pdf](http://www.skcinc.com/pdf/1812.pdf)
8. **Hayes, H., Mitzel, R.:** "AIHce 2011 Presentation: Field Evaluation of Diffusive Samplers for Indoor Air VOC Measurements (.pps)" at [www.skcinc.com/prod/690-101.asp](http://www.skcinc.com/prod/690-101.asp)



**Figure 1. ULTRA Passive Sampler with Blank/Correction Compartment**



**Column:** RTX-5, 30 m, 0.32 ID, 1  $\mu$ m df  
**Injection Temp.:** 250 C  
**Detector Temp.:** 250 C  
**Carrier Gas:** Nitrogen  
**Flow:** N<sub>2</sub> = 2.5 ml/min  
 Air = 450 ml/min  
 H<sub>2</sub> = 45 ml/min  
**GC:** Perkin Elmer Clarus® 500  
**Integrator:** TotalChrom Software  
**Oven Program:**  
*Initial Temp.:* 35 C  
*Initial Hold:* 7 min  
*Ramp 1:* 5 C/min to 100 C hold 0.5 min  
*Ramp 2:* 15 C/min to 175 C hold 0.0 min  
**Total Run Time:** 25.5 min

**Figure 2. Sample Chromatogram – VOCs on Tenax TA Sorbent**