


Low-volume PUF Tubes

For Sampling Semi-volatiles from 1 to 5 L/min



Cat. No. 226-92



Cat. No. 226-124



Low-volume PUF
Tube Holder
Cat. No. 224-29P



Multi-purpose
Calibration Jar
Cat. No. 225-111

- Meet specifications of EPA Method TO-10A, EPA IP-8, and ASTM D4861
- Can be used with personal sample pumps up to 5 L/min
- Precleaned to meet method specifications
- Ready to use
- Suitable for ambient, indoor, or personal air sampling
- Sample up to 24 hours
- Custom PUF/sorbent combinations available



Low-volume PUF Tubes

For Sampling Semi-volatiles from 1 to 5 L/min

Ordering Information

Low-volume PUF Tubes

Methods	Compounds	Sorbent (Amount)	Cat. No.
EPA TO-10A, IP-8	Organochlorine and organophosphorus pesticides, carbamate, pyrethrin, triazine, and urea pesticides	PUF (76 mm)	226-92
ASTM D4861	Organochlorine and organophosphorus pesticides and PCBs	PUF/Tenax/PUF (30 mm/750 mg/30 mm)	226-124
ASTM D4861	Organochlorine and organophosphorus pesticides and PCBs	PUF/Tenax/PUF (30 mm/750 mg/30 mm)	226-143
EPA TO-10A	Organochlorine and organophosphorus pesticides, carbamate, pyrethrin, triazine, and urea pesticides	PUF/XAD-2/PUF (30 mm/1500 mg/30 mm)	226-126
Custom	Pesticides, PCBs, and PAHs	PUF/Glass Fiber filter (76 mm)	
Custom	Pesticides, PCBs, and PAHs	PUF/Glass Fiber filter (76 mm)	

Note: PUF tubes used in the Statements of Work for Superfund sites must be used and analyzed within 14 days of being cleaned.

Accessories	Cat. No.
Low-volume PUF Tube Holder includes fitting with durable protective cover, 3 feet of tubing, and collar clip <i>Do not use an adjustable low flow holder.</i>	224-29P
Multi-purpose Calibration Jar , standard size	225-111
Replacement Parts	
Replacement PUFs for 226-92 , 76-mm length Uncleaned, pk/20	P22692
Cleaned, pk/10	P22692C
Replacement Glass , precleaned and ready to use	P22692G

SKC Limited Warranty and Return Policy

SKC products are subject to the SKC Limited Warranty and Return Policy, which provides SKC's sole liability and the buyer's exclusive remedy. To view the complete SKC Limited Warranty and Return Policy, go to <http://www.skcinc.com/warranty>.

